

Aleavit DAP



COMPOSITION

Dibasic ammonium phosphate $(\text{NH}_4)_2\text{HPO}_4$

GENERAL CHARACTERISTICS

Appearance: whitish powder with faint ammonia odor

Determinazione	U.M.	Intervallo di Conformità	Risultato
Titolo	%	$96,0 < x < 102,0$	Conforme
P ₂ O ₅	%	> 53,5	Conforme
Azoto Ammoniacale	%	> 20,8	Conforme
Sostanze Insolubili	%	< 0,05	Conforme
SO ₃	ppm	< 8130,0	Conforme
Arsenico (As)	ppm	< 2,0	Conforme
Piombo (Pb)	ppm	< 1,0	Conforme
Cloruri (Cl)	ppm	< 10,0	Conforme

APPLICATIONS

ALEAVIT DAP is a fermentation activator designed to create the most favorable environment for the multiplication and growth of yeasts, prolonging their viability and thus avoiding difficult or stunted fermentations. The use of ALEAVIT DAP allows the nitrogen in the must to be supplemented, enabling more regular fermentation even after the first few days.

The use of ALEAVIT DAP helps to avoid the formation of hydrogen sulfide.

RECOMMENDED DOSAGES

10 to 30 g/hl, depending on the APA content of the grapes.

20 to 30 g/hl in refermentation and in fermentation stops.

MODE OF USE

Dissolve Aleavit DAP in a little water and add it in pumping over to the must or wine to be fermented, taking care to properly homogenize the mass. Aerate conveniently as appropriate.

PACKAGES AVAILABLE

25 Kg Bag

STORAGE CONDITIONS

Unopened package: store in a cool (temperature below 25°C), dry and ventilated place.

Open package: close tightly and store as indicated above.

ALEA
Evolution

Via Sandro Pertini, 12
40062 Molinella (BO) - ITALIA
Tel. +39 051 88 7052
+39 346 58 33 882
info@alea-evolution.com
alea-evolution@pec.it
alea-evolution.com
CF - P.IVA 02944251202
Cod. REA BO-479327
Capitale sociale € 40.000,00 i.v.





ALEA
Evolution

Via Sandro Pertini, 12
40062 Molinella (BO) - ITALIA
Tel. +39 051 88 7052
+39 346 58 33 882
info@alea-evolution.com
alea-evolution@pec.it
alea-evolution.com
CF - P.IVA 02944251202
Cod. REA BO-479327
Capitale sociale € 40.000,00 i.v.

